



Environment Testing
America

PFAS: Producing Defensible Data




Taryn McKnight,
PFAS Practice Leader
Eurofins Environment
Testing America


1

Per and Poly?

Perfluorinated = Completely Fluorinated



Polyfluorinated = Incompletely Fluorinated




Environment Testing
America

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

Surfactant Properties

Fluorocarbon "Tail" = Hydrophobic and Oleophobic






Functional Group "Head" = Hydrophilic



Environment Testing
America

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.


Branched & Linear Isomers



Environment Testing
America

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

Replacement Chemicals			
	Legacy Manufacturers		
Original Chemical	PFOA	PFOS	PFOS
Replacement Chemical	HFPO-DA "GenX"	DONA	F-53B



EPA Method 533

“A Method for Short Chain PFAS”



eurofins

Environment Testing
America

533	537.1
Drinking Water	Drinking Water
Branched/Linear Isomers -YES	Branched/Linear Isomers - YES
14 of the same and 11 unique compounds	14 of the same and 4 unique compounds
SPE WAX	SPE SDVB
Hold Time: 28/28 days	Hold Time: 14/28 days
LCMSMS with confirmation ion	LCMSMS - no confirmation ion
Isotope Dilution	Internal standard
Recovery Correction - YES	Recovery Correction – NO
RLs: Not defined	RLs: 2ppt - 40ppt

7

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.


Tissue, Soil, Groundwater, Air?



What methods do we follow for non-potable water, solid and air matrices?


8

Methods Timeline – Past to Present




Environment Testing
America

9
Copyright © 2020, Eurofins Environment Testing America. All rights reserved.



Environment Testing
America

10
Copyright © 2020, Eurofins Environment Testing America. All rights reserved.




11

“PFAS by LCMSMS Compliant with Table B-15 QSM 5.1 or latest version”



12



Analyte Description	NPW & Solids
Perfluorobutanoic acid (PFBA)	EPA Draft Target Analyte List
Perfluoropentanoic acid (PFPeA)	
Perfluorohexanoic acid (PFHxA)	
Perfluoroheptanoic acid (PFHpA)	
Perfluorooctanoic acid (PFOA)	
Perfluorononanoic acid (PFNA)	
Perfluorodecanoic acid (PFDA)	
Perfluoroundecanoic acid (PFUnA)	
Perfluorododecanoic acid (PFDoA)	
Perfluorotridecanoic Acid (PFTriA)	
Perfluorotetradecanoic acid (PFTeA)	
Perfluorobutanesulfonic acid (PFBS)	
Perfluorohexanesulfonic acid (PFHxS)	
Perfluoroheptanesulfonic Acid (PFHpS)	
Perfluorooctanesulfonic acid (PFOS)	
Perfluorodecanesulfonic acid (PFDS)	
Perfluorooctane Sulfonamide (FOSA)	
N-methyl perfluorooctane sulfonamidoacetic acid (NMeFOSAA)	
N-ethyl perfluorooctane sulfonamidoacetic acid (NEtFOSAA)	
Perfluoro-1-pentanesulfonate (PFPeS)	
Perfluoro-1-nonanesulfonate (PFNS)	
6:2FTS	Replacement Chemicals
8:2FTS	
4:2FTS	
DONA	
HFPO-DA (GenX)	
F-53B Major	
F-53B Minor	

Key Elements in Reducing Variability



Environment Testing
America

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

537.1 vs DRAFT 8327 vs 537M

PFAS Method Comparison Table for Aqueous Matrix				
	Features	Method 537.1	EPA DRAFT 8327	User Defined "537 Modified" ¹
	Matrices	Drinking Water	Non-potable aqueous	All aqueous matrices
	Sample size	250 mL	Any allowed (5mL used)	250 mL
1	Analysis	LCMSMS	LCMSMS	LCMSMS
2	Aqueous Extraction	SPE SDVB	DAI	SPE Waters WAX
3	Branched/Linear Isomers	Yes, for available standards	Yes, for available standards	Yes, for available standards
4	Confirmation Ion	No	Yes	Yes
5	Quantitation	Internal standard	External standard	Isotope dilution
	Reporting Limits	(2 ppt - 40 ppt)	(10ppt - 8000ppt)	(2ppt - 20ppt)
	Recovery Correction	No	No	Yes

¹ – Compliant with DoD QSM Table B-15




Environment Testing
America

16
Copyright © 2020, Eurofins Environment Testing America. All rights reserved.



Solid Phase Extraction (SPE)



**SEPARATION
CHEMISTRY**


Branched and Linear Isomers

STANDARD

SAMPLE



Environment Testing
America



19

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.


Secondary Ion Transition - PFOS

STANDARD

SAMPLE



Environment Testing
America



20

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

Isotope Dilution – Labeled Analogues




The Parr Family = Native PFOS


+  =




The Incredible Family = Labeled PFOS



+  =



 = $^{13}\text{C}_{21}$

 eurofins

Environment Testing
America

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

Benefits of Isotope Dilution

What affects the native analyte will equally affect the isotope

Calibration

Most accurate and precise method

Target analytes are quantitated against structurally similar materials, the isotopes themselves

Matrix Mitigation

Expands ability to process a broader range of matrices

Compound Identification

Reduces the potential for false positives

Reduces the potential for error; corrects for retention time shifts

 eurofins


Environment Testing
America

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.


22

EIS response must be > 10:1 S/N

Native = 100 ppt



Native = 100 ppt



Environment Testing
America

23

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

Isotope Dilution vs Internal Standard


PFOA	Spike amount	Isotope Dilution Result	% recovery (ID)	Recovery Corrected Result	Internal Standard Result	% recovery (IS)
Field Sample		10		11	9	
Matrix Spike	34	40	88%		40	91%
Spike Duplicate	34	42	94%		42	97%

MeFOSAA	Spike amount	Isotope Dilution Result	% recovery (ID)	ID Recovery Corrected Result	Internal Standard Result	% recovery (IS)
Field Sample		40		80	40	
Matrix Spike	40	20	50%		20	50%
Spike Duplicate	40	20	50%		20	50%



Environment Testing
America


Copyright © 2020, Eurofins Environment Testing America. All rights reserved.



Environment Testing
America

25

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.



Environment Testing
America

26

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

Managing Artifacts



Field Crew:

personal care, clothing, food, visitors, notebooks, tarps

Sampling Equipment:
verified to be PFAS free?



Sample Collection:


wash hands, wear gloves, don't filter, include field QC



Environment Testing
America

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.

Bottle Types by Collection Method



Environment Testing
America

Robel et al., 2018

28

Copyright © 2020, Eurofins Environment Testing America. All rights reserved.



<https://www.serd-estcp.org/Program-Areas/Environmental-Restoration/Contaminated-Groundwater/Monitoring/ER19-1205>






CA Water Board – POTW PFAS Orders

Treatment System Sampling

- Samples shall be representative and taken before the monitored flow joins or is diluted by any other waste stream, body of water, or substance
- 24 hr sample using a composite sampler is preferred – however, there are challenges and limitations to prevent cross contamination so collecting a grab sample maybe better suited for the facility. An updated sampling guide will be available soon on the Water Board's PFAS website.
- Influent samples shall be taken from locations prior to primary settling
- Effluent samples shall be taken following the completion of all treatment processes and prior to mixing with the receiving waters





Taryn McKnight
PFAS Practice Leader
916.347.6815
taryn.mcknight@eurofinset.com



eurofins

Environment Testing
America

THANK YOU FOR ATTENDING

Copyright © 2020 Eurofins Environment Testing America. All rights reserved.

33