




Total Kjeldahl Nitrogen




Presented by:

Joga Chizer and Nicole Van Aken

Method:

Hach Company Method 10242,
Simplified Spectrophotometric Measurement of
Total Kjeldahl Nitrogen in Water and
Wastewater





BACKGROUND AND FORMS OF NITROGEN

WUDGIWIR QDO#WNQ #RUP XOD

INORGANIC NITROGEN
TOTAL NITROGEN



=

Ammonia
Nitrogen

+

Organic
Nitrogen

KDFK #WNQ #RUP XOD



=

Total
Nitrogen

-


Nitrate
and
Nitrite

DIFFERENT TKN METHODS



Hach Method

Distillation Method



WHY SHOULD I CHOOSE THIS METHOD?


PROS

- 1 Fast
- 2 Uses Equipment Already in the Lab
- 3 Contained (Reagents in Kit, No Outside Chemicals Needed)
- 4 Minimal Training Required

CONS


- 1 Cost
- 2 Dependency on Hach Reagents
- 3 Limited Customization
- 4 Method Deviations Needed for Compliance


METHOD DEVELOPMENT



TIMELINE

FROM TRIAL TO ACCREDITATION






Sample Number	QC template Name	Parameter Name	Corrected Result
Bat-240924-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.959
Bat-240911-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.979
Bat-240815-003-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.04
Bat-240813-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.838
Bat-240813-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.716
Bat-240716-002-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.98
Bat-240704-005-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.919
Bat-240619-012-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.08
Bat-240619-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.08
Bat-240515-004-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.919
Bat-240507-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.08
Bat-240419-006-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.08
Bat-240410-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.02
Bat-240317-002-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.18
Bat-240214-010-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.08
Bat-240214-004-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.06
Bat-240119-005-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.06
Bat-231220-004-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.02
Bat-231212-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.22
Bat-231107-001-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.06
Bat-231019-003-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.08
Bat-230921-003-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.959
Bat-230916-002-ICV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.939
Bat-230824-001-CCV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	0.959
Bat-230605-008-CCV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.33
Bat-230605-007-CCV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.12
Bat-230519-006-CCV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.21
Bat-230506-002-CCV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.2
Bat-230421-007-CCV-01	TKN ICV (1.0 mg/L)	Total Kjeldahl Nitrogen	1.15

CHALLENGES

#1: THE LOW-LEVEL CALIBRATION VERIFICATION (ICV) WAS NOT PASSING CONSISTENTLY

10.2 Calibration Verification

10.2.1 To verify that the instrument is measuring TKN properly, analyze a 1.0 mg/L and 10.0 mg/L NH₃-N standard. Results shall be within 10 percent of the actual value. Perform this calibration verification daily while instrument is in use.




MANUAL CALIBRATION

Read standards at 345nm
and plot in excel

Uses $y=mx+B$

Plugs sample absorbance into
the excel calibration curve



CHALLENGES


#2: METHOD MODIFICATION NEEDED

- A Hach removed Microcap C from the TKN test kit
- ↓
- B 40 CFR part 136 approved method cannot be used
- ↓
- C Alternate Test Procedure (ATP) is not required, but EPA approval is needed

Method listed in 40 CFR

Current Hach Method

Method 10242 rev 1.1	Modifications: Method 10242 rev 1.2
No pH adjustment prior to analysis	Adjust preserved samples to pH 7 with 5 N NaOH
Digestion at 100°C	Digestion at 120°C
Digest for 1 hour	Digest for 30 minutes
Add Micro Cap C after digestion	No Micro Cap C
MDL limit is 0.43 mg/L	MDL limit is 0.11 mg/L





PROCEDURE

PROCEDURE

STEP 1: PH ADJUST YOUR SAMPLES

Samples are preserved during collection

Hach method states "Adjust to pH 7"

Tip: Use your sample as the acid if needed



WHAT QUALITY CONTROL (QC) IS NEEDED?

**Every
Batch
Needs:**

Calibration Verification (AKA Initial Calibration Verification, ICV)

Method Blank (AKA Laboratory Reagent Blank, LRB)

Matrix Spike / Matrix Spike Duplicate

Ongoing Precision and Recovery (AKA Laboratory Control Sample, LCS)


PROCEDURE

STEP 2: QC and SAMPLE PREPARATION


Make dilutions before digesting

Add Solution A and Tablet B

Do NOT invert or mix before digestion



PROCEDURE



STEP 3: DIGEST SAMPLES

20mm cell size needed

Hach only sells heat block with eight (8) vial wells

QC take 6 spots (at least), so multiple digestions need to be done

DIGESTION

Digestion converts all forms of nitrogen to nitrate

The digested sample represents the total nitrogen content




STEP 4: ALIQUOT INTO VIALS

Red Vials

Green Vials

A photograph of a laboratory setting. In the foreground, a white box is prominently displayed. The box has a green and red label. The green part of the label contains the text 'TNT880' and 'Simplified TRN TKA Simplified TKA'. The red part of the label contains the text 'TNT880' and 'Simplified TRN TKA Simplified TKA'. The box is placed on a white surface. In the background, there are various laboratory equipment, including a white container, a black container, and a rack of test tubes.

TKN Kit – TNT880



Solution D



Test Vials

PROCEDURE



STEP 5: ANALYZE SAMPLES

We don't use the TNT auto-read, we read absorbance manually @345nm

Because of that, we prepare a zero standard to zero the machine

(Total Nitrogen) – (Nitrate plus Nitrite)



Fairfield Suisun Sewer District Laboratory
Total Kjeldahl Nitrogen Analysis Worksheet
 Hach Method 10242

Document ID: E-WST-TKN-10242
 Revision: 3.0
 Effective: 08/14/2024

LIMS Batch ID: 250114-001
 Calibration Date:
 Analysis Date:
 Analysis Time:
 Analyst:

Read the red vials first

Green vials are subtracted from red vials

Final concentration calculated by pulling from calibration curve

Digestion		Temp @ Start (°C)		Temp @ End (°C)	
Digestion Batch	Left or Right	Temp @ Start (°C)	Temp @ End (°C)	Temp @ Start (°C)	Temp @ End (°C)
1	Left				
2	Right				
3					
4					

Calibration Curve

Cal Std. Conc. (mg/L)	addition	Digestion Batch	Abs Red Vial	Abs Green Vial	@ 345 nm	Cal Std. Conc.
Cal Blank	1430	1446	1	0.000	0.000	0.000
1.00	1430	1446	1	0.045	-0.004	0.049
3.00	1430	1446	1	0.155	0.000	0.155
7.00	1430	1447	1	0.350	0.000	0.350
10.00	1430	1447	2	0.493	0.000	0.493
15.00	1431	1447	2	0.746	0.000	0.746

Unhide Rows 12-19 to view calibration

LIMS ID	Sample Date & Time	Sample Name	pH after adjust. (6.5 to 7.5)	Digestion Batch	Dilution Factor	Time of Solution addition	Read Time	Abs Red Vial	Abs Green Vial	Absorbance @ 345 nm	Conc. (mg TKN-N/L)	Final Conc. (mg TKN-N/L)	%Rec	RPD
Bat-250114-001-ICV-01		QC-TKN ICV (1.0 mg/L)	N/A	1	1	0940	0955	0.053	0.000	0.053	1.04	1.04	103.96	
Bat-250114-001-ICV-02		QC-TKN ICV (10 mg/L)	N/A	1	1	0940	0956	0.515	0.001	0.514	10.34	10.3	103.42	
Bat-250114-001-MB-01		QC-TKN MB	7.09	1	1	0941	0956	0.003	0.002	0.001	-0.01	-0.01	<	
25JAN02-01-02	1/2/25 7:12	Influent, INF-001	6.76	2	10	1111	1126	0.171	0.003	0.168	3.36	33.6	<	
25JAN02-01-02	1/2/25 7:12	Influent, INF-001MS	6.76	2	1	1111	1126	0.435	0.004	0.431	8.67	8.67	106.13	
25JAN02-01-02	1/2/25 7:12	Influent, INF-001MSD	6.76	2	1	1112	1127	0.441	0.001	0.440	8.85	8.85	109.77	2.07
25JAN13-01-01	1/13/25 7:45	Influent, INF-001	6.67	2	10	1112	1127	0.200	0.004	0.196	3.93	39.3	<	
Bat-250114-001-LCS-01		QC-TKN LCS/OPR	7.31	1	1	0941	0956	0.262	0.002	0.260	5.22	5.22	104.33	
					1								<	



Thank you
Questions ~ Comments ~ Thoughts...



For more information:
JChizer@FairfieldSuisunSewer.ca.gov
NVanAken@FairfieldSuisunSewer.ca.gov

A Special Thank You To:
Jennifer Vo & Raymona Shirmard